



UBL201-L Introductory Biology III - Immunology

Rahul Chavan - 22086

12th September 2023

Quantitation of nitric oxide (NO) production by macrophages in response to different doses of lipopolysaccharide (LPS)

Principle:

Nitric oxide (NO) is a free radical that is produced by the enzymatic conversion of L-arginine to L-citrulline by nitric oxide synthase (NOS). NO is produced by a variety of cells, including macrophages, in response to inflammatory stimuli such as lipopolysaccharide (LPS). NO is a potent vasodilator and has been shown to be involved in the pathogenesis of septic shock. Inflammatory stimuli can enhance NO release via the up-regulation of the inducible form of NOS (iNOS or NOS2) Within inflammatory macrophages. NO is converted into highly Reactive Nitrogen Species (RNS) such as NO_3^- and NO_2^- within infected macrophages to drive bacterial death. NO is also involved in the regulation of apoptosis and the immune response.

NO production can be measured by the Griess reaction, which is based on the ability of NO to react with sulfanilamide and N-(1-naphthyl)ethylenediamine. Sulfanilamide is converted to a diazonium salt by reaction with nitrite (NO_2^-) in acid solution. The diazonium salt is then coupled to N-(1-naphthyl) ethylenediamine leading to the formation of an azo dye (Pink colour), a chromophore that can be measured spectrophotometrically at 540 nm. The amount of NO produced by macrophages in response to LPS can be quantitated by comparing the absorbance of the sample to a standard curve generated using a known concentration of sodium nitrite. Using a standard curve, the concentration of NO produced by the macrophages can be determined from the absorbance of an unknown sample. Therefore, using Griess's reagent we can indirectly estimate levels of NO by quantification of nitrite in the biological samples.




Figure 1: Reaction of nitrite with Griess reagent

Results and discussion:

Absorbance values of serial dilutions of NO₂ at 550 nm:

Table 1: Absorbance values of serial dilutions of NO₂ at 550 nm

Concentration of NO ₂ (μM)	208.3	104.15	52.075	26.0375	13.01875	6.5009	3.2505	unknown
Absorbance Values for Trial 1	2.7295	1.4981	0.811	0.407	0.2853	0.1093	0.0914	0.3657
Absorbance values for Trial 2	3.8812	0.4592	0.4985	0.2243	0.2316	0.115	0.0864	0.3733
Absorbance values for Trial 3	3.8595	0.7958	0.2099	0.1189	0.0979	0.0957	0.0713	0.3616

Absorbance values of the blank at 550 nm:

Table 2: Absorbance values of the blank at 550 nm

Concentration of NO ₂ (μM)	Blank
Absorbance Values for Trial 1	0.0449
Absorbance values for Trial 2	0.0451
Absorbance values for Trial 3	0.00455

In our case the absorbance values for Trials 2 and 3 would give us a parabolic curve and not a straight line. In order to infer for the concentrations, we chose the absorbance values of Trial 1 only, to account for the concentration of the unknown sample.


Concentration vs Absorbance Table:

Table 3: Concentration vs Absorbance Table

Concentration of NO ₂ (μM)	Absorbance
208.3	2.6846
104.15	1.4532
52.075	0.7611
26.0375	0.3621
13.01875	0.2404
6.5009	0.0644
3.2505	0.0465

Here, the new absorbance values are obtained by subtracting the absorbance value of the blank from the original values (values are from Trial 1 only).

Concentration vs Absorbance Graph:



Linear fitting of the plot:




Figure 2: Linear fitting of the plot

Using the plot, we can determine the concentration of the unknown sample:

As we know : $y = mx + b$

where, y = concentration of NO₂ (µM)

x = absorbance (550 nm)

m = slope of the line

b = y-intercept

For the unknown, $x = 0.3208$

from the plot we get, $m = 77.2$ and $b = -2.92$

therefore, $y = 21.84576$

Hence, the nitric oxide concentration of the unknown sample using Griess reagent assay is 21.84576 µM.

Interpretation:

From the plot we can see that there exists a linear relationship between the concentration of NO₂ and the absorbance values. As we know that more is the LPS dosage to the cells, more is nitrite production, we can also establish that more LPS dosage samples will give higher absorbance values. The blank as expected gave negligible absorbance values as it does not contain any nitrite. However, due to pipetting errors, the difference in the absorbance values of duplicates and triplicates is significantly high. Hence, we cannot rely on the values of the duplicates and triplicates. We can only rely on the values of the first trial.

Using the data we can conclude that nitric oxide concentration of the unknown sample using Griess reagent assay is 21.84576 μM i.e., it lies between the 4th and 5th diluted concentrations which are 26.0375 μM and 13.01875 μM respectively, implying the source cells were subjected to a lower concentration of LPS dosage or the cells were not able to produce enough nitric oxide due to lower activation.

Precautions:

- Wear appropriate personal protective equipments, including lab coats, gloves, and safety goggles, to protect against chemical exposure.
- Use sterile techniques and clean glassware to prevent contamination of samples and reagents.
- Calibrate and standardize the spectrophotometer using appropriate standards and controls to ensure accurate measurements.
- Avoid pipetting errors and ensure that the pipettes are calibrated.
- Avoid exposure to light, as Griess reagent is light-sensitive.
- Perform duplicates and triplicates to reduce the error and enhance reproducibility of results.
- Careful disposal of biohazardous waste.